


Digital Imaging and Communications in Medicine DICOM

Pantelis Theodorakis

Bottom- Up Approach



Attribute/Elements (e.g. Patient Name)

TAG	VR	VL	Value
(0010, 0010)	PN	10	Smith^Joe^

- ▶ Tag: (Group #, Element #) to identify the attribute/data element (more than 2000 available)
- ▶ Value Representation (VR): data type to encode the value(s) e.g. PN = Patient Name up to 64 characters
- ▶ Value Length (VL): The length of the Value stored (always even-pad with space if needed)

Byte#	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
Decimal	16	0	16	0	‘P’	‘N’	10	0	S	m	i	t	h	^	J	o	e	(space)
Binary	10	00	10	00	50	4E	0A	00	53	6D	69	74	68	5E	4A	6F	65	20
	g = 0010		e = 0010		VR type		VR length L = 10 = 0x000A		VR value = Smith^Joe (with trailing space)									

Module (e.g. Patient Module)

Attribute	Tag	Type	Attribute Description
Patient Name	(0010,0010)	2	Patient's Full Name
Patient ID	(0010,0020)	1	Primary hospital identification number
Issuer of Patient ID	(0010,0021)	3	Identifier of the Assigning Authority that issued Patient ID
...	

- ▶ Module is a logical group of attributes about a common topic
- ▶ Type: (1) Required (2) May be empty if unknown (3) Optional (C) Conditional

IOD - Information Object Definition (e.g. Enhanced CT Object)

IE	Module	Reference	Usage
Patient	Patient	C.7.1.1	M
...
Equipment	General Equipment	C.7.5.1	M
Image	General Image	C.7.6.1	M
	Contrast/Bolus	C.7.6.4	C-Required if contrast media was used in this image
	CT Image	C.8.2.1	M

- ▶ Information Entity (IE): A group of modules representing a real-world object
- ▶ An IOD is a collection of Modules
- ▶ A Module is a collection of attributes from one IE
- ▶ Usage: (M) Mandatory, (C) Conditional, (U) Optional

DICOM Real World Models

- ▶ A Patient can have (0-n) Visits to a hospital/clinic
- ▶ Physicians order (0-n) Imaging Service Request for each patient
- ▶ Radiology plans (1-n) Requested Procedures which are reported and billed to satisfy each ISR
- ▶ Radiology Information System (RIS) schedules (1-n) Scheduled Procedure Steps for each RP
- ▶ Equipment executes (1-n) Performed Procedure Steps for each SPS


DICOM Information Model

- ▶ A Patient may have many studies
 - ▶ A Study is a group of Series for a given Requested Procedure
 - ▶ A Series may group closely related Images from same PPS, same protocol & same piece of equipment
 - ▶ An Image holds the acquired data
- ▶ Instances are actual data created based on an IOD
- ▶ DICOM uses unique identifier (UIDs) to identify
 - ▶ Specific Instances
 - ▶ Specific SOP Classes
 - ▶ Specific Study/Series
 - ▶ ... and many other things

DICOM Services

- ▶ Print - Printing IODs to a DICOM Printed
- ▶ Storage - Storing Objects, e.g. to Pacs
- ▶ Query/Retrieve - Getting Objects, e.g. from a PACS
- ▶ MWM - Getting Scheduled Patients, e.g. from RIS (Modality Worklist Management)
- ▶ MPPS - Status (Started, Completed) back to RIS (Modality Performed Procedure Step)
- ▶ ... many more


DICOM SOP Class



- ▶ SCP-Service Class Provider
 - ▶ The system that provides the service
- ▶ SCU-Service Class User
 - ▶ The system that uses the service

DICOM Association Negotiation

- ▶ Before two Applications Entities (AE) perform DICOM transaction they must first agree:
 - ▶ Which SOP Class will they use (e.g. MR Image Storage)
 - ▶ Who will be the SCU, who is the SCP
 - ▶ Which is the Transfer Syntax (e.g. Explicit Little Endian)
- ▶ This is the Association Negotiation



Why is it so specific?

- ▶ Complex medical image Data (e.g.)
 - ▶ MR Image
 - ▶ Single frame, 12-16 bit grayscale image
 - ▶ MR acquisition-pulse sequence parameters
 - ▶ 3D patient relative co-ordinate/vector position
 - ▶ X-Ray Angiography Image
 - ▶ Multi-frame, 8-10 bit grayscale image
 - ▶ XA acquisition - radiation/collimation/motion
 - ▶ Dynamic C-arm/table relative positioning

Product DCS

- ▶ DICOM Conformance Statement
 - ▶ Lists the SOPs supported by the product
 - ▶ Describes product implementation details and behaviors
 - ▶ Describes any proprietary attributes

Thank you for your attention!!!